MAP SYSTEM

Direction for use

MAP System®**Intended purpose and Indication**

MAP System is an instrument intended for the placement of the root canal repair material as part of a local obturation procedure.

The local obturation is indicated in the following clinical cases:

- Pulp capping (orthograde application)
- Iatrogenic root perforation (orthograde application)
- Pathological root perforation (root (radicular) resorption) (orthograde application)
- Apexification and apexogenesis of immature teeth (orthograde application)
- Retro-obturation (retrograde application)
- Cervical root resorption (retrograde application)

Properties

The MAP System is composed of a syringe, 5 types of needles available in different sizes and compatible with the respective color-coded plungers. The indications for each type of needle are reported in table ① on page 7.

Intended user

The MAP System is intended to be used only by professional in dentistry field (general dentist practitioners, endodontists).

Patient target groups

The MAP System can be used in individuals of any age, gender or demographic group. There are no demographic restrictions for the use of the MAP System.


Direction for use

1. Clean and sterilize all the components according to the protocol described in section cleaning and sterilization.


2. Insert the same color-coded plastic plunger into the corresponding needle.

3. Press the piston of the MAP Syringe so that the bayonet appears.

4. Lock the plunger inside the bayonet. Then, release the piston of MAP Syringe.


5. Screw the needle onto the syringe. If using a NiTi Memory Shape needle, manually adjust the needle to the desired curvature. Be aware not to bend the needle to close to the base. A space of 1 cm should be left from the base.


6. Prepare the repair material on a glass tab and gather the material into a mound - or into the well supplied in the MAP Kit. Load the end tip of the MAP Needle by tapping several times.

7. Remove excess material by using a gauze.
8. Press the syringe piston to expel the material. Fill the instrument again if necessary.
9. When the material has been fully placed, press the syringe piston 3x right after use and no later than 15 minutes to expel the excess material.

10. If a NiTi Memory Shape needle is used, straighten the needle manually. Unscrew the needle.
11. Press the syringe piston and take the plunger out of the bayonet system.
12. Clean the inside of the needle with a spongy and stiffened-end dental floss like Cleany Floss. If a MAP VPT Needle is used, clean the inside of the needle with the provided MAP Cleaning Brush. They should be entered in the needle in the sense of the needle as shown in picture below.

13. Discard the used plastic plunger and the used dental floss. The MAP Cleaning Brush can be reused.
14. Clean all other parts immediately in order to remove any residual material.

Cleaning and sterilisation

Operation	Operating Mode	First Use	Following Use
Pre-cleaning	Initial treatment at the point of use - Eject the repair material directly after use. The ejection shall be done within 15 minutes after use. - Disassemble the device. Separate the plastic plunger from the needle and discard the plastic plunger. - For needles other than VPT, clean the inside of the needle with a spongy and stiffened-end dental floss such as PD Cleany Floss or EMOFORM TrioFloss. For VPT needles, clean the inside of the needle with the provided MAP cleaning brush. If the dental floss is too thin, several dental flosses can be used at the same time. - Immerse all reusable components in an enzymatic detergent for at least 5 minutes, following the manufacturer's recommended instructions. The detergent should not contain di- or triethanolamines as a corrosion inhibitor (see general recommendations). - Do not use a detergent containing phenol or other substances incompatible with the instruments (see general recommendations). - If visible impurities are observed on instruments, brush the surface using a soft-bristled brush for minimum 20 seconds.		●
Rinsing	Preparation before cleaning Rinse the device with potable water for at least 1 minute.		●
Manual cleaning assisted by an ultrasonic device	Manual cleaning - Place the devices in a kit or container to avoid any contact between instruments. - Prepare an enzymatic detergent solution according to the manufacturer's recommended instructions and place it in a sonication unit. The detergent should not contain di- or triethanolamines as a corrosion inhibitor (see general recommendations). - Completely submerge the device in the cleaning solution and sonicate for a minimum of 15 minutes. - Rinse the device using distilled water for a minimum of 3 minute.	●	●
Rinsing	Drying Allow the device to air dry on a non-linting wipe or with filtered compressed air until no visible traces of moisture remain.	●	●
Inspection	Inspection and maintenance - Visually check the cleanliness of the device. If dirt is observed, repeat the cleaning steps until there are no more signs of dirt. - Inspect devices and sort out those with defects. Discard instruments that have any deformation (bent, twisted), damage (broken, corroded) or defects (loss of color code or marking) affecting the resistance, safety or performance of the instrument.	●	●
Plunger preparation	Pre-sterilisation - MAP Plungers need to be cleaned and sterilized before use. - Clean the MAP plunger with a non-linting wipe soaked in alcohol.	●	
Packaging	Packaging - Do not reassemble the device prior to sterilization. - Place the needle, the plunger and the syringe independently in double sterilization packaging (pouches made of paper and plastic film). - For USA, use FDA-cleared pouches. - Check the validity period of the pouch given by the manufacturer to determine the shelf life. - Physical-chemical indicators may be used on or in the packaging system. - Once packaged, sterilize as soon as possible.	●	●

Operation	Operating Mode				First Use	Following Use																									
Sterilization	Sterilization	- Use a steam sterilizer conforming to EN 13060, EN 285. - Use the sterilization parameters mentioned in table below, which have been validated.																													
		<table><thead><tr><th>Sterilizer Type</th><th>Minimum Temperature</th><th>Minimum exposure Time (min.)</th><th>Minimum Drying Time (min.)</th></tr></thead><tbody><tr><td>Pre-Vacuum</td><td>132°C</td><td>4</td><td>20</td></tr><tr><td>Pre-Vacuum</td><td>134°C</td><td>3</td><td>20</td></tr><tr><td>Pre-Vacuum</td><td>134°C</td><td>5</td><td>20</td></tr><tr><td>Pre-Vacuum</td><td>134°C</td><td>18</td><td>20</td></tr><tr><td>Pre-Vacuum (for plunger)</td><td>132°C</td><td>3</td><td>20</td></tr></tbody></table>				Sterilizer Type	Minimum Temperature	Minimum exposure Time (min.)	Minimum Drying Time (min.)	Pre-Vacuum	132°C	4	20	Pre-Vacuum	134°C	3	20	Pre-Vacuum	134°C	5	20	Pre-Vacuum	134°C	18	20	Pre-Vacuum (for plunger)	132°C	3	20		
		Sterilizer Type	Minimum Temperature	Minimum exposure Time (min.)	Minimum Drying Time (min.)																										
		Pre-Vacuum	132°C	4	20																										
		Pre-Vacuum	134°C	3	20																										
		Pre-Vacuum	134°C	5	20																										
		Pre-Vacuum	134°C	18	20																										
		Pre-Vacuum (for plunger)	132°C	3	20																										
		At the end of the cycle, check the chemical indicator and controls to confirm the efficiency (packaging integrity, no humidity, color change of sterilization indicators, digital records of cycle parameters).																													
Storage	Storage	- Keep devices in sterilization packaging in a dry and clean environment. - Sterility cannot be guaranteed if packaging becomes opened, damaged, or wet. - Check the packaging and the medical devices before using them (packaging integrity, no excessive humidity and validity period).																													

Precautions

To reduce risk of non-functional needles, use only compatible root canal repair materials for the intended purpose.

Disassembly the device before cleaning and sterilisation.

In order to ensure the proper functionality of the device and to avoid pre-mature breaking of MAP NiTi needles, bend the needle leaving a 1 cm space from the base.

Before each use, check appearance of defects such as cracks, deformations (bent, twisted), corrosion, loss of color coding or marking. Those are indications that the device is not able to fulfill the intended use with the required safety level.

Compatible root canal repair materials

MTA, Calcium Hydroxyde Powder and bioceramic cements such as Biodentine™* by Septodont, TotalFill® BC RRM™ Fast Set Putty and TotalFill® BC RRM™ Putty by FKG Dentaire.

*only compatible with blue, white and purple needles.

General recommendations for reprocessing

Use only a disinfecting solution which is approved for its efficacy (VAH/DGHM-listing, CE marking, FDA approval) and in accordance with the instructions of the disinfecting solution manufacturer. It is recommended to use anticorrosion disinfecting and cleaning agents (like: Alconox™ Liquinox™ Critical Cleaning Liquid Detergent).

For your own safety, please wear personal protective equipment (gloves, glasses).

The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments where applicable after sterility.

The device is validated for 50 reprocessing cycles.

The appearance of defects such as cracks, deformations (bent, twisted), corrosion, loss of color coding or marking, are indications that the device is not able to fulfill the intended use with the required safety level.

Plastic and NiTi parts are degraded by Hydrogen Peroxide (H₂O₂) solution. NiTi Instruments are degraded if immersed more than 5 minutes in a solution of NaOCl at more than 5 %.

Do not use acid (pH< 6) or alkaline (pH> 8) solutions.

Don't use any direct flame to straighten the NiTi Needles.

Adverse reactions

None known.

Contraindications

None known.

Warning messages

MAP System needs to be cleaned and sterilized before use.

Clean inside the needle within 15 minutes with a stiff threader dental floss like Cleany Floss or, if MAP VPT Needles are used, with the MAP Cleaning brush. If this is not done, it may be very difficult, if not impossible, to unblock and clean the device.

The MAP Plunger are single-use devices. Do not reuse them in order to avoid risk of cross-contamination between different patients or the use of degraded product. Reprocessing of plastic plunger is not recommended.

Do not use Biodentine with yellow or red needles in order to avoid clogging the needle.

Do not use damaged components.

Storage

Plastic plungers should be kept away from direct light.

Shelf life

The MAP plungers have a shelf life of 5 years. Do not use the product after the expiry date.

Disposal

MAP System components are to be considered contaminated after each use. Non-functional components and plungers should be disposed of according to best practices and local regulations.

In case of serious incident

Any serious incident, leading to temporary or permanent serious deterioration of patient's, user's or other person's state of health, that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Maintenance

In case of breakage of a component, spare parts can be ordered separately.

Kit packaging description**MAP Intro Kit**

- 1 NiTi Memory Shape Needle, No. 0
- 1 NiTi Memory Shape Needle, No. 1
- 1 NiTi Memory Shape Needle, No. 2
- 1 Stainless steel syringe
- 16 Plastic plungers, No. 0
- 16 Plastic plungers, No. 1
- 16 Plastic plungers, No. 2
- 1 Stainless steel well

Sold in procedure pack with Cleany Floss

MAP One Kit

- 1 NiTi Memory Shape Needle, No. 1
- 1 Stainless steel syringe
- 16 Plastic plungers, No. 1
- 1 Stainless steel well

Sold in procedure pack with Cleany Floss

MAP Surgery Kit

- 1 Surgery Needle, No. 0, right
- 1 Surgery Needle, No. 0, left
- 1 Surgery Needle, No. 1, right
- 1 Surgery Needle, No. 1, left
- 1 NiTi Memory Shape Needle, No. 0
- 1 NiTi Memory Shape Needle, No. 1
- 1 Stainless steel syringe
- 24 Plastic plungers, No. 0
- 24 Plastic plungers, No. 1
- 1 Stainless steel well
- 1 Aluminium box 185x145x30 mm

Sold in procedure pack with Cleany Floss

MAP VPT Kit

- 1 White Needle, No. 1
- 1 Stainless steel syringe
- 16 White plungers, No. 1
- 1 Stainless steel well

Sold in procedure pack with MAP Cleaning brush

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