Deutsch  
English  
Español  
Français  
Italiano

Scheren

Scissors

Tijeras

Ciseaux

Forbici

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Instructions for Use

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Mode d'emploi

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**CE 0297**

|                       |           |
|-----------------------|-----------|
| <b>Deutsch .....</b>  | <b>7</b>  |
| <b>English .....</b>  | <b>25</b> |
| <b>Español .....</b>  | <b>47</b> |
| <b>Français .....</b> | <b>68</b> |
| <b>Italiano .....</b> | <b>89</b> |

**Symbolerklärung**  
**Symbol Explanation**  
**Explicación de los símbolos**  
**Explications des symboles**  
**Spiegazione dei simboli**

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## Table of Contents

|          |                                                                              |           |
|----------|------------------------------------------------------------------------------|-----------|
| <b>1</b> | <b>General Information .....</b>                                             | <b>27</b> |
| 1.1      | Manufacturer .....                                                           | 27        |
| 1.2      | Hotline .....                                                                | 27        |
| 1.3      | Event Reporting Duty .....                                                   | 27        |
| 1.4      | Notes on this Document .....                                                 | 27        |
| 1.5      | Terms & Acronyms .....                                                       | 28        |
| 1.6      | Validity of this Document .....                                              | 28        |
| 1.7      | Accompanying Documents .....                                                 | 28        |
| <b>2</b> | <b>Scope of Delivery .....</b>                                               | <b>28</b> |
| 2.1      | Inspecting the Delivery for Completeness and Correctness .....               | 28        |
| <b>3</b> | <b>Intended Use .....</b>                                                    | <b>29</b> |
| 3.1      | Intended Purpose .....                                                       | 29        |
| 3.2      | Indications .....                                                            | 29        |
| 3.3      | Contraindications .....                                                      | 29        |
| 3.4      | Target Patient Groups .....                                                  | 29        |
| 3.5      | Users .....                                                                  | 29        |
| 3.6      | Ambient Conditions During Application .....                                  | 29        |
| 3.7      | Warnings .....                                                               | 29        |
| <b>4</b> | <b>Application .....</b>                                                     | <b>30</b> |
| 4.1      | Description of the Components .....                                          | 30        |
| 4.1.1    | Structure, Operating Principle and Performance Specifications .....          | 30        |
| 4.1.2    | Materials Used .....                                                         | 32        |
| 4.1.3    | Variants / Combinations .....                                                | 32        |
| 4.1.4    | Compatible Products and Accessories .....                                    | 32        |
| 4.2      | Before Initial Use .....                                                     | 33        |
| 4.3      | Application .....                                                            | 33        |
| 4.4      | After Conclusion of the Surgery .....                                        | 33        |
| <b>5</b> | <b>Cleaning, Disinfection and Sterilization .....</b>                        | <b>34</b> |
| 5.1      | Limitations and Restrictions to Cleaning, Disinfecting and Sterilizing ..... | 34        |
| 5.2      | Pre-Treatment at the Application Location prior to Cleaning .....            | 35        |
| 5.3      | Preparation before Cleaning .....                                            | 35        |
| 5.3.1    | Preparation before Cleaning for Instruments and Storage Systems .....        | 36        |
| 5.4      | Cleaning and Disinfection .....                                              | 36        |
| 5.4.1    | Manual Cleaning and Disinfection .....                                       | 36        |
| 5.4.2    | Machine Cleaning and Disinfection .....                                      | 37        |
| 5.5      | Inspection, Functional Check, Maintenance .....                              | 39        |

|          |                                               |           |
|----------|-----------------------------------------------|-----------|
| 5.5.1    | Inspection and Functional Check .....         | 39        |
| 5.5.2    | Maintenance .....                             | 41        |
| 5.6      | Packaging.....                                | 41        |
| 5.7      | Sterilization .....                           | 41        |
| 5.8      | Storage and Transport.....                    | 42        |
| <b>6</b> | <b>Maintenance.....</b>                       | <b>43</b> |
| 6.1      | General Notices .....                         | 43        |
| 6.2      | Accessories & Spare Parts .....               | 43        |
| <b>7</b> | <b>Ecological Information / Disposal.....</b> | <b>44</b> |
| 7.1      | Packaging.....                                | 44        |
| 7.2      | Disposal .....                                | 44        |
| 7.3      | National Regulations .....                    | 44        |

## 1 General Information

### 1.1 Manufacturer

We thank you for buying one of our products.

This product bears the CE mark. This means that it meets the basic requirements on the safety and the performance of medical devices set forth by the respective applicable European regulations.

We are the manufacturer of this product:

**KLS Martin SE & Co. KG**

A company of the KLS Martin Group

KLS Martin Platz 1 · D-78532 Tuttlingen · Germany

Phone +49 7461 706-0 · Fax +49 7461 706-193

[info@klsmartin.com](mailto:info@klsmartin.com) [www.klsmartin.com](http://www.klsmartin.com)

### 1.2 Hotline

Should you have any questions on how to handle the device or product or use it for clinical applications, please do not hesitate to contact the Product Management:

Phone: +49 7461 706-0

E-mail: [info@klsmartin.com](mailto:info@klsmartin.com)

### **NOTICE**

Each packaging and most products are marked with a specific batch code (LOT) and a catalog number (REF). Please always provide LOT and REF in case of a complaint.

### 1.3 Event Reporting Duty

All serious adverse events occurring in connection with the product must be reported to KLS Martin and the respective authorities without delay.

### 1.4 Notes on this Document

### **WARNING**

**Possible danger to the life of patients, users and others when these Instructions for Use are disregarded!**

Read Instructions for Use completely and follow them carefully. In particular, be sure to note all cautionary and warning notices.

This document refers to persons of all genders alike. Reference to various genders is avoided purely for reasons of improved readability.

The electronic version of these Instructions for Use can be requested at <https://www.klsmartin.com/en/services/instructions-for-use/>.

## 1.5 Terms & Acronyms

| Acronym        | Description                                                                         |
|----------------|-------------------------------------------------------------------------------------|
| A <sub>0</sub> | Standard for the killing of microorganisms in procedures with humid heat, ISO 15883 |
| AKI            | Instrument processing working group                                                 |
| CJD            | Creutzfeldt-Jakob disease                                                           |
| OR             | Operating room                                                                      |
| WD             | Washer-disinfector                                                                  |

## 1.6 Validity of this Document

This document applies to the product family of reusable Scissors.

## 1.7 Accompanying Documents

The specific information for every reference number for the section “Cleaning, disinfection, and sterilization” is also applicable with these instructions for use. It is available at:

[www.klsmartin.com/processing](http://www.klsmartin.com/processing)

## 2 Scope of Delivery

The scope of delivery depends on the product selected and may include various items allocated to the article group described under section 1.6 "Validity of this Document".

### 2.1 Inspecting the Delivery for Completeness and Correctness

Immediately upon receipt, the goods must be checked for completeness and damages.

Notice of any such damage must be given immediately.

After delivery of the product, check the original packaging and the packaging seal for damages.

## 3 Intended Use

### 3.1 Intended Purpose

The scissors, saws and drills are reusable surgically non-invasive instruments that are not used in connection to an active product and whose function is in the transient cutting, piercing, severing and drilling of dressing materials or other medical auxiliary materials during a treatment.

### 3.2 Indications

Surgical treatment of pain or diseases.

Because of the intended use of the products, all surgical interventions and regions of the human body are indicated.

### 3.3 Contraindications

There are no known contraindications when the product is used safely.

### 3.4 Target Patient Groups

There are no restrictions regarding the patient target group.

### 3.5 Users

The products must only be used by qualified surgical staff.

The cleaning, disinfecting and sterilizing is to be performed by trained and qualified staff of the medical device preparation unit.

### 3.6 Ambient Conditions During Application

Use only in ORs under surgical conditions or within specified zones.

### 3.7 Warnings

**Possible danger to the life of others by shipping contaminated products!**

In case of return shipments, send only clean and disinfected products in sterile packaging.

---


## 4 Application

### 4.1 Description of the Components

#### 4.1.1 Structure, Operating Principle and Performance Specifications


Detailed information can be found in the special brochure produced by KLS Martin.

##### Surgical scissors





- |                         |               |
|-------------------------|---------------|
| 1 Cutting edges (bevel) | 3 Branches    |
| 2 End, lubricate here   | 4 Ring handle |

##### Dissection scissors

- |                         |           |
|-------------------------|-----------|
| 1 Cutting edges (bevel) | 3 Handles |
| 2 End, lubricate here   | 4 Spring  |

Surgical scissors, dissection scissors, and microscissors are available in the following designs:

| Figure                                                                                                      | Description                                  |
|-------------------------------------------------------------------------------------------------------------|----------------------------------------------|
| <br><b>sutureCut</b>       | Suture material scissors (ligature scissors) |
| <b>finoLine</b>                                                                                             | Slim model                                   |
|                            | straight                                     |
|                            | curved                                       |
|                            | blunt/blunt                                  |
|                          | sharp/sharp                                  |
|                          | sharp/blunt                                  |
|                          | Stainless steel                              |
| <br><b>TC GOLD</b>       | TC hard metal in the working end             |
| <br><b>SuperCut</b>      | Scissor cut or knife cut                     |
| <br><b>TCC BlackLine</b> | TC hard metal with knife cut                 |

| Figure                                                                            | Description                              |
|-----------------------------------------------------------------------------------|------------------------------------------|
|  | TC hard metal with knife cut or wave cut |
| <b>A</b>                                                                          | for Adventitia                           |
| <b>D</b>                                                                          | for Dissectia                            |

#### 4.1.2 Materials Used

| Designation              | International name                                           |
|--------------------------|--------------------------------------------------------------|
| Surgical stainless steel | 1.4021 / 1.4024 / 1.4034 / 1.4117 / 1.4301 / 1.4305 / 1.4310 |
| Titanium                 | 3.7165                                                       |
| Plastics                 | PEEK/PTFE                                                    |
| Coatings                 | gold-coated/anodized/TiAlN (SolidBlack)/AlCrN + CBC          |
| Further additives        | Tungsten carbide                                             |

#### 4.1.3 Variants / Combinations

The products are available in various variants. Please consult the special brochure produced by KLS Martin for this.

The products can be found in the following catalogs at the minimum:

| Designation                | REF          |
|----------------------------|--------------|
| General, C35               | 90-100-48-07 |
| Oral maxillofacial surgery | 90-112-48-05 |
| Heart surgery              | 90-132-48-04 |
| Dental                     | 90-138-48-08 |
| Plastic surgery            | 90-272-48-04 |
| Neurosurgery               | 90-424-48-04 |

#### 4.1.4 Compatible Products and Accessories

No specific system accessories are available for this product.

## 4.2 Before Initial Use

Remove any protective caps or protective coverings required for transport before using the instrument for the first time.

We recommend washing all new surgical steel instruments separately before introducing them into the instrument cycle for positive impact on the formation of the passive layer. Factory-new surgical steel instruments in particular are more susceptible to surface changes as a result of the impact of products that are already part of the process.

After this initial step, the products can be incorporated in the usual reprocessing routine. In order to do this, follow the instructions from chapter 5 “Cleaning, Disinfection and Sterilization”.

## 4.3 Application

### Opening the scissors

Ring handle: Pull the rings apart and the working part will open.

Spring handle: When the spring is relaxed, the working part is already open.

### Cutting

Ring handle: Push the rings together and the working ends will glide over one another and the scissors will cut.

Spring handle: Push the handle surfaces together and the working ends will glide over one another and the scissors will cut.

## 4.4 After Conclusion of the Surgery

---

### **NOTICE**

Even improper depositing can cause damage to the products. Care should therefore be taken when handling the product.

---

Deposit the products properly and transfer to reprocessing in accordance with chapter 5 “Cleaning, Disinfection and Sterilization”.

## 5 Cleaning, Disinfection and Sterilization

### **Non-sterile handling can result in danger to the life of the patient!**

Products supplied in non-sterile condition must be cleaned, disinfected and sterilized before the first and before any further application.

- In case of patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible CJD variants, reprocessing of the products must be performed in compliance with applicable national regulations.
- The operator / preparer is responsible for the cleaning, disinfecting and sterilizing of the used instruments. It is essential that national regulations, including restrictions, be observed.

The preparation process must be validated before the medical device is used. The operator / preparer is responsible for this.

The specific information as per the reference number of the following section is available at:

[www.klsmartin.com/processing](http://www.klsmartin.com/processing)

### **5.1 Limitations and Restrictions to Cleaning, Disinfecting and Sterilizing**

The alloying of steels ("stainless") used in the manufacturing of instruments forms special passive protection layers.

These steels have only limited resistance to chloride ions and aggressive media and liquids!

In addition to KLS Martin's efforts in choosing and carefully processing the correct materials, surgical instruments must be continually and professionally maintained and correctly processed at the hands of the user. Furthermore, the medical device must be subjected to the previously validated proper preparation.

Frequent reprocessing has little effect on the surgical instruments.

The product lifetime of a surgical instrument is primarily determined by wear and possible tear during use.

This is detected during the functional and visual inspection to be performed during the individual reconditioning cycles. If necessary, appropriate corrective measures must be initiated. This functional and visual inspection is the prerequisite for the medical device's release for the next reconditioning, for the storage of the medical device in the validated sterile barrier system and for the application of the medical device after performed sterilization.

A general statement with respect to a maximum of cycles is therefore not possible.

Advice on functional and visual check procedures can be found in section 5.5, "Inspection, Functional Check, Maintenance"

## 5.2 Pre-Treatment at the Application Location prior to Cleaning

Conduct the cleaning and disinfecting of contaminated instruments as soon as possible after use.

The following generally applies:

The sooner the internal lumina are cleaned, the better the result will be. Should it not be possible to comply with this urgency due to the length of the application or due to organizational aspects, it is the user's sole responsibility to specify and validate measures in order to prevent the drying of the contaminants or to ensure the success of the cleaning despite drying.

Generally, we recommend return transport in dry condition inside a closed transport container without adding fluids, cleaning detergents or disinfectants.

## 5.3 Preparation before Cleaning

---

### CAUTION

#### **Possible danger of injury to the patient caused by residues on the medical device!**

Products with rough surfaces, threads, sharp edges, narrow gaps or similar features can capture particles from cloths and brushes.

- Use only soft, clean and lint-free cloths for drying.
- The brushes for ducts must be a little bit bigger than the inner diameter of the duct and the shaft of the brush must be at least as long as the duct.

---

A machine-based process (WD) must be used for cleaning and disinfecting. Whether a product is validated for a manual procedure can be reviewed at [www.klsmartin.com/processing](http://www.klsmartin.com/processing). Due to the substantially lower efficiency and reproducibility, a manual process – even when using an ultrasonic bath – should only be applied when a machine-based process is not available.

When choosing a cleaning detergent, the compatibility with the material as well as the suitability for and the efficiency in the cleaning of medical devices must be considered.

The manufacturer instructions for the cleaning or combined cleaning and disinfecting agent regarding concentrations, temperature, application time as well as rinsing specifications must be observed. Use only freshly prepared solutions, only sterile or low-germ water (max. 10 germs/ml) and only low endotoxin load (max. 0.25 endotoxin units/ml) water (e.g. purified water / highly purified water). For drying, use only a soft, clean and lint-free cloth and/or filtered air.

Place the medical devices on carriers suitable for cleaning, e.g. sieve trays or tray baskets. Observe the WD manufacturers' operating and loading instructions.

### 5.3.1 Preparation before Cleaning for Instruments and Storage Systems

The pre-treatment process was validated according to the following procedure:

1. Open or dismantle products to the extent possible.
2. Rinse product under running water for at least 1 min. (temperature < 35°C/95°F).
3. Place the dismantled or opened products into an adequately sized pre-cleaning bath<sup>1</sup> (in a ultrasonic bath) for the prescribed application time, ensuring that the products are fully submerged. Ensure that products do not touch each other. Support the pre-cleaning by fully brushing all surfaces (at the beginning of the application time). Move movable parts back and forth at least three times during pre-cleaning.  
Rinse all lumina of the products at least three times at the beginning and at the end of the application time (aids and minimum volumes depend on the cavity to be rinsed).
4. Where necessary, activate the ultrasonic bath for another minimum application time period (at least 5 min.).
5. After that, remove the products from the pre-cleaning bath and rinse thoroughly with water at least three times (at least 1 min.). Move movable parts back and forth at least three times during rinsing.  
Rinse all lumina of the products at least three times (aids and minimum volumes depend on the cavity to be rinsed).

## 5.4 Cleaning and Disinfection

### 5.4.1 Manual Cleaning and Disinfection

When choosing the used cleaning agents and disinfectants, it must be ensured,

- that they are generally suitable for the cleaning and disinfecting of invasive medical devices made from metals and plastics,
- that the cleaning agent is suitable for ultrasound cleaning (no foaming),
- that a disinfectant with tested efficacy (e.g. VAH/DGHM or FDA/EPA approval / clearance / registration or CE marking) is used and that it is compatible with the used cleaning agent and
- that the used chemicals are compatible with the product. When choosing the cleaning agents and disinfectants, it must be ensured that they do **not** contain the following ingredients:
  - organic, mineral and oxidizing acids (minimum permissible pH value 5.5)
  - organic solvents (e.g. alcohols, ether, ketones, benzines)
  - oxidants (e.g. hydrogen peroxide)
  - halogens (chloride, iodine, bromine)
  - aromatic/halogenated hydrocarbons

---

<sup>1</sup> When cleaning- and disinfecting agents are used, e.g. for occupational safety reasons, they should be aldehyde-free (otherwise fixation of blood and contaminations) with proven efficiency (e.g. VAH/DGHM or FDA/EPA approval / clearance / registration or CE marking). The cleaning and disinfecting agent should be suitable for the disinfection of the products. The sole purpose of the disinfectant used in pre-treatment is personal safety. It cannot replace the process step of disinfection after cleaning!

For products made of aluminum or other materials that are alkali-sensitive, use only alkalis/strong alkalis/neutral/enzymatic agents with a maximum pH-value of 8.5 or alkaline cleaning agents with a maximum pH-value of 11.

Combined cleaning and disinfecting agents should not be used, if possible. Combined cleaning and disinfecting agents can only be used in the case of very minor contaminations (no visible soiling).

When manually cleaning with a present danger of injury or infection, additional occupational safety measures (e.g. protective clothing, protective goggles, gloves, ambient air filtering) must be applied in accordance with national requirements.

The proof of the general suitability of the products for an effective manual cleaning and disinfection was provided by an independent, accredited and recognized test lab (§ 15 (5) MPG [German Medical Devices Act]) using the pre-cleaning and cleaning agents Cidezyme/Enzol and the disinfectant Cidex (Johnson & Johnson GmbH, Norderstedt). In this process, product-specific procedures were observed.

The specific information about the procedures for each reference number are available at:

[www.klsmartin.com/processing](http://www.klsmartin.com/processing)

#### 5.4.2 Machine Cleaning and Disinfection

When choosing a WD, it must be ensured that,

- the WD efficiency is proven (e.g. DGHM or FDA approval / clearance / registration or CE marking according to DIN EN ISO 15883),
- if possible, a tested program for thermal disinfection ( $A_0$ -value > 3000 or – for older devices – at least 5 min. at 90 °C/194 °F) is used (with chemical disinfection there is a risk of disinfectant residues on the products),
- the used program is suitable for the products and has enough rinsing cycles.

Detergent residues can effectively be prevented with any of the following measures:

- a) performing at least three diluting steps after cleaning or neutralization or
- b) conductance control of the WD
- using only sterile or low-germ (max. 10 germs/ml) as well as low endotoxin load (max. 0.25 endotoxin units/ml) water (e.g. purified water / highly purified water)<sup>2</sup>,
- ensuring that the air used for drying is filtered (oil-free, low germ- and particle load) and
- ensuring that the WD is regularly serviced, inspected and calibrated.

---

<sup>2</sup> If a lower water quality should be deemed sufficient in light of national recommendations (e.g. in Germany KRINKO/RKI/BfArM recommendation for preparation), the responsibility falls solely on the preparer.

When choosing the cleaning agent system it must be ensured,

- that it is generally suitable for the cleaning of invasive medical devices made from metals and plastics,
- that, when no thermal disinfection is applied, an additional suitable disinfectant with tested efficacy (e.g. VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) is used and that it is compatible with the used cleaning agent and
- that the used chemicals are compatible with the product.

When choosing the cleaning agents and disinfectants, it must be ensured that they do **not** contain the following ingredients:

- organic, mineral and oxidizing acids (minimum permissible pH value 5.5)
- organic solvents (e.g. alcohols, ether, ketones, benzines)
- oxidants (e.g. hydrogen peroxide)
- halogens (chloride, iodine, bromine)
- aromatic/halogenated hydrocarbons

For products made of aluminum or other materials that are alkali-sensitive, use only alkalis/strong alkalis/neutral/enzymatic agents with a maximum pH-value of 8.5 or alkaline cleaning agents with a maximum pH-value of 11.

The manufacturer instructions for the cleaning agent regarding concentrations, temperature, application time as well as rinsing specifications must be observed.

**Procedure:**

1. Insert dismantled or opened products into the WD. Ensure that products do not touch each other.  
Enable active flushing by connecting to the rinse connector of the WD.
2. Start program.
3. Disconnect the products (if required) after program end and remove from the WD.
4. Inspect and package products as soon as possible after removal.

The following cleaning process has been validated:

| Phase | Step                 | Temperature          | Time                 | Water quality       | Additional Information                                                            |
|-------|----------------------|----------------------|----------------------|---------------------|-----------------------------------------------------------------------------------|
| 1     | Pre-cleaning         | cold (not tempered)  | 1 min.               | Drinking water      | –                                                                                 |
| 2     | Cleaning             | 50 °C                | 5 min.               | Drinking water      | Minimum concentration as specified by cleaning agent manufacturer                 |
| 3     | Neutralization       | –                    | –                    | –                   | Neutralization not carried out according to detergent manufacturer's instructions |
| 4     | Final rinsing        | cold (not tempered)  | 1 min.               | Demineralized water | –                                                                                 |
| 5     | Thermal disinfection | 90 °C                | 5 min.               | Demineralized water | A <sub>0</sub> value ≥ 3000                                                       |
| 6     | Drying               | according to program | according to program | –                   | Drying was not applied during validation                                          |

The proof of the general suitability of the products for an effective machine cleaning and disinfection was provided by an independent, accredited and recognized test lab (§ 15 (5) MPG [German Medical Devices Act]) using the WD G 7836 CD (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the pre-cleaning and cleaning agent neodisher MediClean (Dr. Weigert GmbH & Co. KG, Hamburg). The above-mentioned procedure was observed in doing so.

Metallic coated surfaces (e.g. SolidBlack) and all titanium components may show material changes after only a short while if alkaline cleaning agents or acidic neutralizers are used. These are usually grayish or brownish minor discolorations of the coated surface or a fading of the anodized titanium. According to current research, this does not affect the functionality. Hydrogen peroxide must not be used on SolidBlack coatings. A high degree of material preservation is achieved by using pH-neutral treatment agents.

## 5.5 Inspection, Functional Check, Maintenance

### 5.5.1 Inspection and Functional Check

#### **NOTICE**

Further use is confirmed by the successful inspection of the product. The inspection release and the packaging into a sterile barrier system release the product for its next use.

- The products must be macroscopically clean, i.e. free from visible contamination, after each cleaning process.
- At the end of work, check the instruments, and in particular their connections, for fractures, cracks, deformation, damage and proper function.
- DIN 96298-3 (Medical Instruments – Terms, measuring methods and tests, Part 3: Testing) can offer guidance during functional testing.


- Worn, corroded, deformed, porous or otherwise damaged products must be replaced. Alternatively, appropriate measures can be initiated (e.g. surface treatment, repair), see also the recommendations by the Arbeitskreis Instrumenten-Aufbereitung (AKI, Instrument Processing Working Group):

<http://www.a-k-i.org>

The following evaluation criteria for the functional check and visual inspection of surgical instruments are based on many years of experience. Follow the stated action recommendations. These are performed at the KLS Martin Service.

The following evaluation criteria are an aid to the detection of defects on the product or similar products. A comprehensive evaluation matrix is available from KLS Martin under REF 90-460-01-04.

Products that show one of the defects listed below must not be used.

| Functional limitation                |                                                                                     |                                           |
|--------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------|
| Movement stiff/smooth                |   | Repair: Align                             |
| Blunt/ distorted                     |  | Repair: Re-sharpening of the working part |
| TC worn / chipped                    |  | Exchange                                  |
| Corrosion                            |                                                                                     |                                           |
| Friction corrosion/contact corrosion |  | Repair: Surface treatment                 |
| Pitting corrosion                    |  | Exchange                                  |
| Surface corrosion                    |  | Repair: Surface treatment                 |
| Black surface                        |                                                                                     |                                           |
|                                      |  | Repair: Surface treatment                 |

| Fracture               |                                                                                    |                           |
|------------------------|------------------------------------------------------------------------------------|---------------------------|
| In the pivot           |   | Exchange                  |
| General surface damage |                                                                                    |                           |
| Scratches              |   | Repair: Surface treatment |
| Silicate stains        |   | Repair: Surface treatment |
| Water stains           |   | Repair: Surface treatment |
| Gold worn              |  | Repair: Surface treatment |

### 5.5.2 Maintenance

“Care” means treating the instruments with instrument oil or milk (white oil-in-water emulsion). The joints of the instruments must be treated with steam-sterilizable, paraffin-oil-based care agents after cleaning but before sterilization.

LubriPen®, REF 55-997-01-04 by KLS Martin for the maintenance of surgical instruments.

### 5.6 Packaging

Approved sterilization packaging (e.g. conforming to EN 868, ISO 11607) must be used for sterilization, subsequent transportation and storage.

### 5.7 Sterilization

Sterilization is to be done applying the following steam sterilization method. Other sterilization methods are not permissible.

### Steam sterilization:

- fractionated vacuum process (with adequate product drying<sup>3</sup>)
- Steam sterilizer according to DIN EN 13060/DIN EN 285 or ANSI AAMI ST79 (for USA: FDA-Clearance)
- must be validated according to DIN EN ISO 17665 (valid IQ/OQ (commissioning) and product-specific performance evaluation (PQ))
- Maximum sterilization temperature 134 °C (273 °F, plus tolerance according to DIN EN ISO 17665)

| Steam sterilization with fractionated pre-vacuum |                                 |                                 |                        |
|--------------------------------------------------|---------------------------------|---------------------------------|------------------------|
|                                                  | Procedure 1                     | Procedure 2                     | Procedure 3            |
| Pre-vacuum cycles                                | min. 3 ×                        | min. 3 ×                        | min. 3 ×               |
| Medium                                           | saturated water steam           | saturated water steam           | saturated water steam  |
| Sterilization time                               | 5 min. (or longer) <sup>4</sup> | 5 min. (or longer) <sup>4</sup> | 4 min. (or longer)     |
| Sterilization temperature                        | 134 °C, plus tolerance          | 132 °C, plus tolerance          | 132 °C, plus tolerance |
| Drying time                                      | 20 min <sup>3</sup>             | 20 min <sup>3</sup>             | 20 min <sup>3</sup>    |

The proof of the general suitability of the products for an effective steam sterilization was provided by an independent, accredited and recognized test lab (§ 15 (5) MPG [German Medical Devices Act]) using the steam sterilizer HST 6x6x6 (Zirbus technology GmbH, Bad Grund), applying the fractionated vacuum process and using a commonly available paraffinic-white-oil-based, additive-free instrument oil (greasing of joints and friction surfaces). In doing so, the typical conditions at hospitals and medical offices as well as the above-mentioned procedure were taken into consideration.

Do not use gravitation sterilization, no flash sterilization, no hot air sterilization, no radiation sterilization, no formaldehyde or ethylene dioxide sterilization and also no plasma sterilization.

## 5.8 Storage and Transport

- Store and transport cleanly.
- Avoid severe temperature fluctuations.
- Protect them against mechanical damage.

In case of return shipments, send only clean and disinfected instruments in sterile packaging.

<sup>3</sup> The actual required drying time depends directly on the parameters, which are the sole responsibility of the user (loading configuration and density, sterilizer condition, ...), and must therefore be determined by the user. Drying times should, nevertheless, not be less than 20 min.

<sup>4</sup> or extended sterilization time (e.g. 18 min.) for the prion activation according to national standards

## 6 Maintenance

### 6.1 General Notices

The product may be repaired only by KLS Martin or a qualified person or firm expressly authorized by KLS Martin to perform such work.

Modification of the product may lead to unforeseen risks and is therefore not permitted.

If the repair is carried out by a person or firm specially authorized by KLS Martin, the operator of the product is required to obtain from the repairer a certificate with details about the nature and scope of the repair work done. This certificate must show the date of the repair and the details of the person or firm carrying out the work and must be signed.

In all cases where a party other than KLS Martin performed the work, repaired products must be additionally marked with the repairer's ID label.

Improper interventions or alterations performed by third parties during the period of limitation shall void any and all warranty claims. Unauthorized actions performed on the product are never allowed and shall void any claims of liability against KLS Martin.

For additional information about the maintenance services of the KLS Martin Group, please visit

[www.klsmartin.com/de/services/instrumentenmanagement/](http://www.klsmartin.com/de/services/instrumentenmanagement/)

or contact us directly at

[marmanagement@klsmartin.com](mailto:marmanagement@klsmartin.com)

### 6.2 Accessories & Spare Parts

See the special brochure produced by KLS Martin.

The products can be found in the following catalogs at the minimum:

| Designation     | REF          |
|-----------------|--------------|
| General, C35    | 90-100-48-07 |
| Heart surgery   | 90-132-48-04 |
| Dental          | 90-138-48-08 |
| Plastic surgery | 90-272-48-04 |
| Neurosurgery    | 90-424-48-04 |

## 7 Ecological Information / Disposal

### 7.1 Packaging

KLS Martin will take back the full packaging upon request. Whenever possible, parts of the packaging will be reused.

When disposing of or recycling the packaging, please take applicable national regulations into account.

### 7.2 Disposal

In designing the product, we tried to avoid using composite materials wherever possible. This allows a high degree of recycling. We therefore offer you to take the product back for proper disposal and recycling.

### 7.3 National Regulations

The national regulations and disposal provisions must be observed for all disposal measures.

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